


HOW TO VIEW QRS IMAGES FROM USB DRIVE

1. Flip tab and plug in USB drive (ensure chip side is in correctly)




2. Open "File Explorer" and click on USB drive




3. Double click on patient folder




4. Double click on "Viewer"



5. Double click on "StartIntelViewerCD"



6. This should bring up the following screen, double click on the patient to view images



Dicom Dir

1 Studies

Location	Patient Name	Patient Birth Date ▾	Patient
	TEST, VISAGE MR		10.202

7. Most studies should load within 30 seconds, however studies with a lot of images (eg. CT, MRI Brain) may take up to 2 minutes to load

